lumiraDx® INR Quality Controls

For Professional Use Only
SPEC-30172 Rev5
ART-00041 Rev5

The LumiraDx INR Quality Controls (hereafter referred to as Quality Controls) are optional liquid quality controls to be used with the LumiraDx Instrument (hereafter referred to as the Instrument) and the LumiraDx INR Test (hereafter referred to as INR Test).

Read these instructions thoroughly before using the Quality Controls.

Inspect the Quality Controls packaging and contents for damage before use. Report any damage to LumiraDx Customer Services and do not use the kit if any damage is observed to the contents. Do not perform Quality Controls testing if the Quality Controls material has not been fully reconstituted.

To ensure that you are using the Instrument, the INR Test and the Quality Controls correctly, read the appropriate Platform User Manual, INR Test Product Insert and in this entire pack insert. In addition, please watch the LumiraDx Platform Training video available at kc.lumiraDx.com. The Quality Controls are intended for professional use only.

The LumiraDx Logo, LumiraDx are trademarks of the LumiraDx Group.

Intended Use

The LumiraDx INR Quality Controls are intended for liquid quality control testing performed on the LumiraDx Instrument when used with the LumiraDx INR Test Strip. The Quality Controls provide users with assurance that the device is performing within specification.

Summary and Explanation of the Test

The LumiraDx INR Quality Controls are an optional quality control for the Instrument when used with the LumiraDx INR Test. The Quality Control material is a device intended for medical purposes for use in a test system to estimate test precision and to detect systematic analytical deviations that may arise from reagent or analytical instrument variation and may be used for proficiency testing.

Reagents

Each Quality Control kit contains lyophilized human plasma with varied levels of coagulation factors, stabilizers, preservatives and a calcium chloride Reconstitution Solution. The Quality Control ranges are assigned with the LumiraDx INR Test Strip. The LumiraDx INR Test Strip performance is traceable to the WHO international Reference standard rIF/16.

Warnings and Precautions

For *in vitro* diagnostic use.

- The plasma control contains human source material that was tested and found nonreactive for HIV antibody, Hepatitis B Surface Antigen and Anti-HCV at the donor stage. This product, as with all human based specimens, should be treated as potentially infectious and handled with proper laboratory safety procedures to minimize the risk of transmission of infectious disease.
- All components of this kit can be discarded as Biohazard waste according to the local guidelines.
- Refer to the product safety data sheet for risk and safety phrases and disposal information.
- The product safety data sheet is available for users upon request. Requirements of the appropriate licensing or accrediting body should be incorporated into your quality control program.
- Exercise the normal precautions required for handling all laboratory reagents.

Storage and stability

- Store controls between 2°C and 8°C (36 - 46°F). DO NOT FREEZE.
- Allow the Quality Controls material to come up to room temperature (15 - 30°C; 59 - 86°F) for at least 10 minutes (but not above two (2) hours) before preparing its reconstitution.
- Unopened, lyophilized controls that are stored between 2°C and 8°C (36 - 46°F) can be used until the expiration date.
- Controls are stable for 30 minutes between 15°C and 30°C (59 - 86°F) after adding the Reconstitution Solution. Do not refrigerate after reconstitution

Carton Contents

3 vials Level 1 Quality Control lyophilized plasma (normal INR levels) (Blue Cap)

3 vials Level 2 Quality Control lyophilized plasma (therapeutic INR levels) (Red Cap)

6 vials of Reconstitution Solution (White Cap)

6 large bulb pipettes

6 small bulb pipettes

INR Quality Controls Pack Insert

- Materials required but not provided with the Control Carton:**
 - LumiraDx Instrument
 - LumiraDx INR Test Strips
 - LumiraDx Connect- If connectivity required (refer to LumiraDx Connect User Manual
 - Quality Control Ranges Product Insert as included in the INR Test Strip Carton

1. Getting ready to test

You will need the LumiraDx Instrument and the following supplies:

- LumiraDx INR Test Strip(s)
- LumiraDx INR Quality Controls Level 1 or Level 2
- 1 Vial of Reconstitution Solution
- Large bulb pipette
- Small bulb pipette

2. Preparing the Quality Control

- Remove the Quality Control vial and Reconstitution solution vial from the pack. Allow the vials to come up to room temperature (15 - 30°C; 59 - 86°F) for 10 minutes (but not above two (2) hours).
- Remove the screw-cap and rubber stopper from the Quality Control vial.
- Open the Reconstitution Solution vial and using the large bulb pipette transfer the entire contents to the selected Quality Control vial. Dispose of pipette.

- Close the Quality Control vial with the screw-cap and swirl it gently ten (10) times using a circular motion. DO NOT shake.
- Let the vial sit undisturbed for at least 4 minutes to ensure that the dried material dissolves completely.

- Before use gently swirl the Quality Control vial again another ten (10) times.

- Make sure the Quality Control is fully reconstituted before use (check visually).

- Use the reconstituted Quality Control solution within 30 minutes between 15°C and 30°C (59 - 86°F) from the time the Reconstitution Solution is added. Do not refrigerate after reconstitution.

Handling the LumiraDx INR Test Strips

To ensure that you are using the INR Test and the Instrument correctly, read the appropriate INR Test Strip Product Insert and Platform User Manual.

Procedure/Performing a Test

Consult the LumiraDx Platform User Manual for instructions on how to analyze a Quality Control sample. Open the foil pouch of the INR Test Strip just before use and insert the test Strip into the LumiraDx Instrument. The instrument will indicate when ready for the sample to be applied.

- Draw up the reconstituted **Quality Control solution into the small bulb pipette**.
- Apply the Quality Control solution to the already inserted INR Test Strip.** Hold the small bulb pipette over the Scan Area/ Application Area of the Test Strip and dispense the Quality Control solution. The LumiraDx Instrument will indicate sample is detected with an audible alert (if the Instrument sounds are enabled). The screen of the LumiraDx Instrument will request the user to close the door. Dispose of pipette.
- Do not open the door while the test is in progress.** The touchscreen will indicate test progress.
- The Result** will appear on the touchscreen within 3 minutes of applying the Quality Control solution and starting the test.

- NOTE: If you need to repeat a test, use a new Test Strip.
- Do not open the door while the test is in progress.** The touchscreen will indicate test progress.
- The Result** will appear on the touchscreen within 3 minutes of applying the Quality Control solution and starting the test.

- NOTE: If you need to repeat a test, use a new Test Strip.
- Expected results**

The Instrument displays the result. Consult the Quality Control Ranges Product Insert in the INR Test Strip carton to confirm the result is within the acceptable range. The result is automatically saved in the memory of the Instrument. The system is working properly and all handling has been done correctly when the test results obtained are within the acceptable control range as written in the Quality Control Ranges Product Insert.
- Unacceptable results**

An out-of-range result is indicated if the result is not within the acceptable range as written in the Quality Control Ranges Product Insert. To resolve out-of-range results or error messages, check for the following:

 - Controls may be expired or were stored improperly.
 - You may not have allowed the Quality Control material to come up to room temperature (15 - 30°C; 59 - 86°F) before use.
 - The control may not have been used within 30 minutes after adding the Reconstitution Solution.
 - The Quality Control may not have been fully reconstituted before use.
 - Make sure you run the test within 15 minutes of removing the Test Strip from its foil pouch.
 - You may not be doing the test correctly. Repeat the control test, using a new Test Strip. Carefully follow the instructions in the Platform User Manual.
 - If you follow these guidelines and still require assistance, call LumiraDx Customer Services at 08000 58647239
 - Refer to the How-to-Use training video on kc.lumiraDx.com

- Performance characteristics**

Quality Control precision was determined in measurement system analysis for the Quality Controls with the LumiraDx INR Test. The results were generated over multiple days, by multiple operators and instruments. The results are shown in INR.

Quality Control Level 1	Mean	SD	%CV	N
	1.25	0.081	6.5	235
Quality Control Level 2	Mean	SD	%CV	N
	4.78	0.310	6.5	234

Limitations

This product is designed as a liquid control for quality control performed exclusively by the LumiraDx INR Platform. The quality control result is subjected to the limitations of the test Platform. Deviations may indicate potential problems with one or more components in the test Platform. The LumiraDx Instrument and INR Test Strips have on-board controls to detect errors and prevent false INR results when analysis is performed. Therefore deviations observed when testing with the LumiraDx INR Quality Controls would not invalidate previous results obtained from LumiraDx INR Tests.

LumiraDx Customer Service
If you have questions call the **LumiraDx Customer Services at 08000 58647239**.

Any adverse results experienced with the use of this product, and/or quality problems should also be reported to LumiraDx Customer Services at the above number or by email at customerservices@lumiraDx.com.

The customer international free phone number will be free from most landlines in countries where the LumiraDx Platform is commercially available, but could be classified as an international call from a mobile and therefore chargeable. The charge will depend on the country you are calling from.

For Return Policy

If there is a problem with the LumiraDx INR Quality Controls you may be asked to return them. Before returning product please obtain a return authorization number from LumiraDx Customer Services. This return authorization number must be on the shipping carton for return. For ordinary returns following purchase, please contact Customer Services for terms and conditions.



Warranty - Limited Warranty

LumiraDx Quality Controls – As per shelf life.
For the applicable warranty period, LumiraDx warrants that each product shall be (i) of good quality and free of material defects, (ii) function in accordance with the material specifications referenced in the pack insert, and (iii) approved by the proper governmental agencies required for the sale of products for their intended use (the "limited warranty"). If the product fails to meet the requirements of the limited warranty, then as the customer's sole remedy, LumiraDx shall either repair or replace, at LumiraDx's discretion, the LumiraDx INR Quality Controls. Except for the limited warranty stated in this section, LumiraDx disclaims any and all warranties, express or implied, including but not limited to, any warranty of merchantability, fitness for a particular purpose and non-impingement regarding the product. LumiraDx's maximum liability with any customer claim shall not exceed the net product price paid by the customer. Neither party shall be liable to the other party for special, incidental or consequential damages, including, without limitation, loss of business, profits, data or revenue, even if a party receives notice in advance that these kinds of damages might result.

The Limited Warranty above shall not apply if the customer has subjected the LumiraDx INR Test Strips and Controls to physical abuse, misuse, abnormal use, use inconsistent with the LumiraDx Platform User Manual, INR Test Product Insert or INR Quality Controls Pack Insert, fraud, tampering, unusual physical stress, negligence or accidents. Any warranty claim by Customer pursuant to the Limited Warranty shall be made in writing within the applicable Limited Warranty period.

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Symbols Glossary

	Temperature limitation
	Manufacturer
	In vitro diagnostic medical device
	Catalogue Number
	Lot Number
	Use-by Date – indicates the date after which the unopened IVD/Quality Control material can not be used.
	*CE Mark*: This product fulfills the requirements of the European Directive 98/79/EC on in vitro diagnostic medical devices.
	Indicates a control material that is intended to verify the performance characteristics of the LumiraDx Instrument. An indication of the level of control is shown in the box, e.g. 1 or 2.
	Indicates that potential biological risks are associated with the Quality Control material.
	Refer to instructions for use.
	Authorized representative in the European Community.

Manufacturer information

LumiraDx UK Ltd, Dumyat Business Park, Alloa FK10 2PB, UK

This pack insert and its content is copyright of LumiraDx Group Limited, 2019 - © All rights reserved, worldwide. Content should be used for use of the LumiraDx products and services only and in line with instructions provided. You may not, except with our express written permission, distribute or commercially exploit the content. Nor may you transmit it or store it in any other form of electronic retrieval system.

 LumiraDx AB, Västra Vågen 5A, 16961 Solna, Sweden

 CE mark applies to LumiraDx Instrument, Test Strips, Quality Controls and Connect Hub only.

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Il presente foglietto illustrativo e il relativo contenuto sono protetti dal copyright di LumiraDx Group Limited, 2019 - © Tutti i diritti riservati in tutti i Paesi al mondo.