



CRP Test Specifications

For *in vitro* diagnostic use.

Intended use

The LumiraDx CRP Test is intended for use with the LumiraDx Platform. It is an automated *in vitro* diagnostic test for near-patient testing, for use by healthcare professionals, for the quantitative determination of C-Reactive Protein in human whole blood (capillary finger stick and venous) and plasma samples. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammation disorders, and associated disease. It is intended for use in patients 18 years of age or older.

Test description

CRP is one of the most prominent acute phase proteins (APP), a protein whose serum concentration increases or decreases during acute or chronic inflammatory conditions. As such it has become a universal biomarker of infection and inflammation for a number of diseases and pathophysiological conditions.

The LumiraDx CRP Test is an easy to use, fast microfluidic immunoassay designed to rapidly quantify CRP levels in whole blood and plasma.

The Test delivers rapid, quantitative results from a single direct fingerstick in 4 minutes at the point of care.

The LumiraDx CRP Test is traceable to the ERM® - DA474/IFCC certified global reference materials.

Method comparison

The method comparison was performed using 3 Test Strip lots with plasma samples (Lithium Heparin) sourced from patients presenting with symptoms of respiratory illness, inflammation, or injury, at hospital Emergency Departments (ED), acute medical units or out-patient clinics. A comparison of 320 CRP measurements with the LumiraDx CRP Test to the RCRP Flex® assay on the Siemens Dimension® Xpand® Plus Integrated Chemistry System analysed by Passing Bablok regression yielded the following statistics: Slope = 1.054, Intercept = -0.66, $r = 0.99$.

Matrix Equivalency

A study was carried out with 40 subjects presenting with symptoms of respiratory illness, inflammation or injury. Samples of capillary fingerstick blood (direct application and transfer tube) and paired whole blood (Lithium Heparin) and plasma (Lithium Heparin) samples were collected and tested and data was analysed by Passing Bablok regression. The results showed equivalency across all matrix types.

Built-in Quality Controls

The LumiraDx Platform is integrated with several built-in quality control checks to ensure that the Instrument and Test Strips are functioning correctly for every test run. These checks include:

- Electrical component operation, heater operation, battery charge state, mechanical actuators, sensors and optical system performance
- Test Strip positioning
- Test Strip expiry
- Monitoring of Test Strip performance and microfluidics controls during test runtime
- The CRP Test contains an Onboard Quality Control (OBC) assay
- Sufficient sample volume
- Hematocrit determination on the test strip

CRP Quality Controls

LumiraDx Multi Quality Controls come in two levels and are available from LumiraDx to complete Quality Control assessment of the Instrument and CRP Test Strips.

Precision

A precision study was carried out in heparinised venous plasma on a protocol based on CLSI EP5-A3. The study was carried out at 3 concentrations of CRP, each was tested in 1 run of 5 replicates per day, for five days across 3 sites. The results of the precision study are summarised in the below table:

CRP concentration (mg/L)	Within Day Precision (%CV)	Between Day Precision (% CV)	Between Site Precision (% CV)	Total precision (% CV)	n
11.9	5.1	1.1	1.7	5.5	75
19.7	4.2	4.3	2.6	6.6	75
123.6	4.9	0.0	3.4	6.0	75

CRP Test Strip Specifications

Displayed results	CRP (mg/L to 1 dp)
Storage temperature	2–30°C (36–86°F)
Operating temperature	15–30°C (59–86°F)
Measuring range	5.0 - 250.0 mg/L Each laboratory should investigate the transferability of the expected values to its own patient population and, if necessary, determine its own reference ranges.
Hook effect	No hook effect is observed with the LumiraDx CRP Test at CRP concentrations up to 987 mg/L.
Minimum sample size	20 µL
Sample type	<ul style="list-style-type: none">• Direct fingerstick• Venous Whole Blood (Lithium Heparin)• Plasma (Lithium Heparin)
Time to result	4 minutes

For more information visit lumiradx.com or contact the LumiraDx Customer Services by email: customerservices@lumiradx.com


This specification sheet and its content are copyright of LumiraDx Group Limited © 2022. All rights reserved, worldwide. Content should be used for the use of the LumiraDx products only and in line with instructions provided. You may not, except with our express written permission, distribute or commercially exploit the content. Nor may you transmit it or store it in any other form of electronic retrieval system other than for the purpose of use of the LumiraDx Instrument or LumiraDx Test Strips. Information provided is subject to change without notice.

Product is not available in all countries and regions. Please check with your local LumiraDx sales representative or distributor for availability in specific markets. Not currently available in the USA.

Manufactured by:
LumiraDx UK Ltd,
Dumyat Business Park,
Alloa, FK10 2PB, UK.
Registration Number:
09206123

**Authorized Representative
in the European Union:**
LumiraDx AB
Västra Vägen 5A
16961 Solna,
Sweden

lumiradx.com



S-COM-ART-00329 R1