

IHE Integration Statement			
Vendor	Product Name	Version	Date
LumiraDx	Consult	5.2	3/1/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Node	None
Consistent Time	Time Client	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Demographic Query HL7 V3	Patient Demographics Supplier	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Cross-reference Manager	None
Cross-Community Access	Initiating Gateway	None
Cross-Community Access	Responding Gateway	None
Cross Community Patient Discovery	Initiating Gateway	None
Cross Community Patient Discovery	Responding Gateway	None
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Registry	None
Cross-Enterprise Document Sharing	Document Registry	Patient Identity Feed (HL7 V3)
Cross-Enterprise Document Sharing	Document Repository	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)